SetUp OpenEphys

How to launch Open-ephys :

- if you want to look at the former data already recorded (perceval or others..) :
 - go to x2goclient session on hubel (where is all the data we need (perceval)
 - do : `/home/programs/open-ephys/open-ephys`
 - it will open the app like that (note that is it important the numbers of plugin on the right, some are vital for the process of the data) :

- add a file reader and we are looking for the dat in the doc :
hubel/data-131/perceval/**Raw/**
 - then click until you find in the file you want to analyze a structure.oebin
- if you want to record and have the detection live :
 - open-ephys on your computer, be sure to have the latest version


How to get the right signal chain

The signal chain is where you can put your plugins (things on the left of the screen that you can slide on the bottom of the screen).

Please follow the instructions and use the same parameters. You can see here the result :


and (la suite) :

1. There is 2 choice for the first plugin :
 - a. **“File Reader”** : if you want to read data already recorded (reminder : if so you have to be on the open-ephys of hubel to get the data from hubel). You selected your file (**taking the Raw then the structure.oebin of the session** (list of possible session : **sleepm tachment sleep tachment sleeps**) corresponding to an hour of three day of recording)

- b. OTHERWISE you will use the plugin **“Acquisition Board”** to get the data
2. **“Channel map”** : **ordonner les channels** :
- Get the .xml done from the excel of Antoine that precises which channel corresponds to which one in “reality” and gives the packages of channels that are in the NR, the PFC and the Hippocampus
 - remark : to have the right order, you should do recording of x minutes then study it to do the excel that Antine is doing to have the translation between which channel number is where. And this before doing the real detection (pre-experience!)
 - and then you use this conversion to write the channel mapping in the .json here is an example of how the plugin channel mapping is reading it : “channel_mapping_Rat004_karadoc” at
hubel/media/data-103/Nathan/Avalanche_project/aval_data/
 - use one program python ? to get the file .json of the channel mapping
 - be aware of which group begins in the list (what are the 32 first channels ? by default should be the ones of the NR)
 - be aware of which channel is not used
 - be aware of which channel correspond to the accelerometers components
 - you can now import this file in the channel map plugin
 - remark : open-ephys starts the counting at 1 BUT in .py or other doc it starts at 0 (so there is a +1 on the open-ephys interface at each number or position of channels)**
 - number = “the name of the channel”**
 - position = “the position regarding to other channels”**
3. if you are in case 1.b : **“Record Node”** to get the data for post analysis (Antoine knows how to do it)
- IMPORTANT REMARK** : where to put the Record Node ? Before or after the channel mapping ? Do you want your raw datas to be in order in advance
4. **“Bandpass Filter”** : filter the channels that you are interested in. ATTENTION : **do not filter the channels corresponding to the accelerometer component !!** 300Hz to 6000Hz.
5. **“Spike Detector”** : Here are the parameters STD 4, take all the channels from the NR


#	Name	Type	Channels	Thresholds	Waveform
1	Electrode 1	SE	[55]	STD: 4.0	FULL
2	Electrode 2	SE	[50]	STD: 4.0	FULL
3	Electrode 3	SE	[63]	STD: 4.0	FULL
4	Electrode 4	SE	[49]	STD: 4.0	FULL
5	Electrode 5	SE	[52]	STD: 4.0	FULL
6	Electrode 7	SE	[53]	STD: 4.0	FULL
7	Electrode 8	SE	[54]	STD: 4.0	FULL
8	Electrode 9	SE	[51]	STD: 4.0	FULL
9	Electrode 10	SE	[56]	STD: 4.0	FULL
10	Electrode 11	SE	[57]	STD: 4.0	FULL
11	Electrode 12	SE	[58]	STD: 4.0	FULL
12	Electrode 13	SE	[59]	STD: 4.0	FULL
13	Electrode 14	SE	[60]	STD: 4.0	FULL
14	Electrode 15	SE	[61]	STD: 4.0	FULL
15	Electrode 16	SE	[64]	STD: 4.0	FULL
16	Electrode 16 (1)	SE	[62]	STD: 4.0	FULL


6. **“Mean Spike Rate”** : Select all the channels that you have analyzed in the Spike Detector (hence from the Nucleus Reunien).

- a. Time constant : 100ms. Output is where (which channel) the mean spike rate will be written.
 - b. ATTENTION : this is overwriting on other data, be sure that the output channel is not one already used.
 - i. maybe adding a new channel is possible ?
 - ii. BUT if we put mean spike rate **after** record node, does the output of meanspike rate (channel that its position is written in the plugin mean spike rate in "Output") overwrite the channel where it is? For instance i am putting my mean spike rate on the channel at the position 33 (that we suppose was carrying other data before doing that), is this channel overwritten by those new values and we will recording those new values instead of the formers ones ? Or as the record node is **before** the mean spike rate in the signal chain there is no issue with that and we will keep recording the original values of this channel ?
7. "TTL TOGGLE PANEL" : if it is not here : CTRL+P and install it and its sister TTL DISPLAY PANEL they should appear together in your plugin list.
- a. this is your interface as an user of the detector of InfraSlow Rhythm. When you click on a number it created a TTL event. While you have not turned it off (by clicking again), the event is still considered as "True" (it is important).
 - b. You can click on **nb 10** to start the integration of the thresholds, and it will after 20 sec of integration start to detect avalanches and Infralow.
 - i. **This is a key point in the detection : the moment where you decide to start the detecting program. Those first 20 seconds after your click are vital because determine the calibration of the detection.**
 - ii. **Thus when you should click** : when those conditions are all fulfilled at the same time :
 1. the rat is asleep and is not moving
 2. When you look at the mean Spike rate and it seems that it is currently in an infralowrhythm (succession of ON/OFF)
 - iii. After having clicked on, **be sure to click off the button in less than 20 seconds**. Other wise the program will restart on and on the integration after the 20 seconds.
 - iv. **If it seems that the integration went badly and thus the detection is bad too : click again on 10 and it will restart the integration and the detection.**
 - c. You can click on **nb 15** to stop Infralow (put the event Infralow from True to False) and demand to rebegin the detection to this point.
8. "Python Processor" : again it is not here CTRL+P etc.. Select the program : **NathPie_detector_user_mode.py** from
 hubel/media/data-103/Nathan/Avalanche_project/aval_code/ (you can copy it locally)
- a. be sure to read the python Processor Section
 - b. Normally you should get your Python Environment by selecting the folder "miniconda3" at
 Users>Acquisition\miniconda3\envs\oe-python-plugin

9. **“LFP Viewer”** : use the split screen to have on the top the visual of the channels and on the bottom the visual of the only channel corresponding to mean spike rate (ie the output channel of this plugin)
10. **“TTL Display Panel”** : to know which event is True or False at the very moment.
 - a. Something good to know : You can use on the bottom screen the Overlay Events to see filters that show if currently we are in avalanche (green nb 8), InfraSlow (violet nb5) and Wake/REM (orange nb 2)

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
How to modify/get access to the Python Processor code:

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

Setup Opto System

How to connect O-E and Arduino


- charge your open-ephys locally (on your computer)
 - add the signal chain corresponding (cf SetUp Open-ephys)
 - /\ open-ephys locally won't work if you try to read files because they are on hubel and mounting large amount of data is too much and graphically very bad.
 - read files of already recorded session on hubel :
/home/programs/open-ephys/open-ephys
- add to the end of the signal chain :
 - the ARDUINO OUTPUT plugin :

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REMARQUE : if you want to measure with oscilloscope : put it on REF OUT of the cyclops to see the result of the voltage you are creating.

SetUp Codes

How to interpret the DATA : Matlab :

- 2 codes exist on Matlab with 2 functionalities at the location :
hubel/media/data-103/Nathan/Avalanche_project/aval_code/
- first code : **Sessions_analyse_UsAvals_Nath.m**
 - to get the DATA from an animal, one day and analyse on session (of a phases like for instance "sleepm")
 - to look at mean firing rate/ IS period / avalanches for this session
 - determine the appropriate start time to use in the Open Ephys File Reader
- second code : **Comparison_OE_Matlab_detection.m**
 - to get the DATA from an animal, one day and analyse on session (of a phases like for instance "sleepm")
 - to get the timestamps of your detection from the open-ephys programme
NathPie_detector_user_mode
 - Temporally align Open Ephys detections with IS and avalanche intervals obtained from Pietro's MATLAB algorithm
 - Sometimes you need **to manually remove the last interval** in **IS_OE_intervals** and **wakeREM_OE_intervals**. In some cases, the Open Ephys output contains only the start timestamp (first column) without a corresponding end timestamp. This results in an invalid interval where the computed end time equals the **start_sec_decay** offset and becomes smaller than the start time, which is not physically possible. Such incomplete intervals must be discarded before running the comparison analysis.
 - Compute statistical metrics to quantitatively compare Open Ephys detection with MATLAB (Pietro) detection

NathPie_detector : usage

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- Where to find it : **hubel/media/data-103/Nathan/Avalanche_project/aval_code/**
- Its name : **NathPie_detector_user_mode.py**
 - a. copy it on your local machine (your computer), not not modify the official version
- If it is not clear in english : **NathPie_detector_user_mode_inFrench.py**
(commentaires originaux)

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SETUP :

1. cf Setup Open-ephys for the signal chain to setup in open-ephys gui
2. open the python programme (on vscode for instance)
3. check the parameters at the section "modifiable parameters"
 - a. see what is next for more precision
 - b. modification depending of the rat / channels (electrode and accelerometer) / computer
4. **not forget to save your changes on your computer**
5. then on open-ephys : fetch your program with the plugin python processor
 - a. if there is a problem the open-ephys will tell you
6. run the acquisition with the recording

HOW IT WORKS :

Processor operates through several successive stages:

1. First, **data acquisition** and moving average computation are performed during 20 seconds (self.time_limit_integration, it is a parameter free to be changed if necessary). Spike rate values and acceleration signals are extracted from the selected channels. A moving average is then computed over the last seconds defined by **time_value_moy** for spike rate and **time_value_acceleration** for acceleration. These smoothed signals are used for all subsequent detections in order to reduce high-frequency fluctuations and improve stability.
2. **User control** is implemented via the TTL toggle panel (cf open-ephys plugin). One TTL line, **ligne_threshold_integration** (default value 9 on python , corresponding to the **10** on the TTL PANEL), is dedicated to triggering threshold recalculation. When the user activates this line, a new threshold integration phase is initiated and all parameters of the last integration disappear to make place for new thresholds. Another TTL line, **ligne_stop_IS_now** (default value 14 on python , corresponding to the **15** on the TTL PANEL), allows the user to immediately stop IS detection. When this line is activated, all ongoing detections and related states are reset.
 - a. VITAL PART : where the user should press the button 10 :
 - i. **Those first 20 seconds after your click are vital because determine the calibration of the detection**
 - ii. **Thus when you should click** : when those conditions are all fulfilled at the same time :

1. the rat is asleep and is not moving
 2. When you look at the mean Spike rate and it seems that it is currently in an infraslowrhythm (succession of ON/OFF)
 - iii. After having clicked on, be sure to click off the button in less than 20 seconds. Other wise the program will restart on and on the integration after the 20 seconds.
 - iv. If it seems that the integration went badly and thus the detection is bad too : click again on 10 and it will restart the integration and the detection.
3. **Threshold computation** occurs after accumulating values during `time_limit_integration` seconds. The spike rate threshold (used after to detect whether or not we are in ON or OFF avalanche) is defined as the percentile specified by `percentage`. The acceleration threshold is defined as the percentile specified by `percentage_accel` (again to define on the accelerometer data whether or not the rat is awake or else different than nREM). In addition, a baseline acceleration level at rest is computed and used later for Wake/REM state detection.
4. **Two things :**
- a. **Wake/REM detection** is triggered when the moving average acceleration exceeds the acceleration threshold or when the spike rate becomes excessively high. In this case, the system switches to a Wake/REM state. Wake/REM periods are recorded, and any ongoing IS or ON detection is reset during these periods to prevent false IS classification. If we are not in the state we can enter in the Avalanche detection and InfraSlow detection
 - b. **Avalanche detection** is based on transitions between OFF and ON states. The transition from OFF to ON is triggered when the moving average spike rate exceeds the spike rate threshold. At this moment, the beginning of the avalanche is recorded in `start_time_ON`, and the previous OFF interval is closed by storing `end_time_OFF`. The transition from ON to OFF occurs when the moving average spike rate falls below the threshold. The end of the avalanche is recorded in `end_time_ON`. If the avalanche duration exceeds `max_time_ON`, it is considered invalid and IS detection is reset. If the duration is within acceptable limits, the corresponding OFF interval duration is computed and stored in `duration_OFF`. But if an ON period is too small we consider it as noise detection and we still are in OFF period
5. **InfraSlow (IS) detection** is based on the temporal structure of OFF and ON intervals. An IS period is detected when at least two OFF intervals (thus 2 ON) are longer than `min_time_OFF`, and none of the OFF intervals exceeds `max_time_OFF` (the selection of max and min ON period is already done in 4.). When these conditions are satisfied, the IS state is updated and the corresponding TTL events are sent. Wake/REM state flags are also updated accordingly.

THE PARAMETERS :

1. The parameter `time_limit_integration` defines the duration in seconds during which spike rate and acceleration values are accumulated to compute detection thresholds. For example, a value of 20.0 means that thresholds are estimated from the first 20 seconds of data (or from the moment threshold recalculation is triggered). A longer integration time produces a more statistically stable threshold but delays detection readiness. A shorter integration time makes the system more reactive but potentially less robust to noise.
2. The parameter `percentage` defines the percentile used to compute the spike rate threshold. For example, a value of 40 means that the spike rate threshold corresponds to the 40th percentile of the spike rate distribution during the integration window. Lower values make the system more permissive, increasing ON detections, whereas higher values make detection more selective.
 - a. For now i have tested this parameter to be optimized for best detection (cf “Reglage du modèle”) at 40
3. The parameter `percentage_accel` defines the percentile used to compute the acceleration threshold. For example, a value of 70 means that the acceleration threshold corresponds to the 70th percentile of the acceleration distribution during integration. Increasing this value makes Wake/REM detection more conservative, while decreasing it increases sensitivity to movement.
 - a. For now i have tested this parameter to be optimized for best detection (cf “Reglage du modèle”) at 70
4. The parameter `max_time_ON` defines the maximum allowed duration, in seconds, for an avalanche (ON state). For example, 14 seconds means that any ON period longer than 14 seconds will be considered physiologically invalid and ignored. When this limit is exceeded, IS detection is reset. This parameter prevents pathological or sustained high-activity periods from being misclassified as valid avalanches.
 - a. For now i have tested this parameter to be optimized for best detection (cf “Reglage du modèle”)
5. The parameter `min_time_ON` defines the minimum duration required for an avalanche to be considered valid. For example, 0.3 seconds means that very brief threshold crossings shorter than 300 ms will be ignored. This prevents noise-driven fluctuations from being classified as true ON events.
 - a. For now i have tested this parameter to be optimized for best detection (cf “Reglage du modèle”)
6. The parameter `max_time_OFF` defines the maximum allowed duration for an OFF interval. For example, 16 seconds means that if an OFF period exceeds 16 seconds, it breaks the temporal structure required for InfraSlow detection. This parameter ensures that IS periods are composed of relatively regular OFF intervals and not separated by excessively long silent phases.
 - a. For now i have tested this parameter to be optimized for best detection (cf “Reglage du modèle”)
7. The parameter `min_time_OFF` defines the minimum duration required for an OFF interval to be counted toward IS detection. For example, 0.5 seconds means that very short OFF gaps are ignored. This avoids counting micro-interruptions in activity as structurally meaningful OFF periods.

- a. For now i have tested this parameter to be optimized for best detection (cf "Reglage du modèle")
8. The parameter `time_value_moy` defines the time window, in seconds, used to compute the moving average of the spike rate. For example, 8 seconds means that each decision of whether or not there is a too high intensity of spike is based on the mean spike rate over the last 8 seconds. A larger window increases stability but reduces temporal precision. A smaller window increases responsiveness but may introduce variability.
 - a. For now i have tested this parameter to be optimized for best detection (cf "Reglage du modèle")
9. The parameter `time_value_acceleration` defines the time window used to compute the moving average of acceleration. For example, 8 seconds means that movement detection is based on smoothed acceleration over 8 seconds. As with spike rate smoothing, larger values increase stability while smaller values increase sensitivity to rapid changes.
 - a. For now i have tested this parameter to be optimized for best detection (cf "Reglage du modèle")
10. The parameter `multiplier_tresh` is a multiplicative factor applied to the spike rate threshold. For example, a value of 3 means that the computed percentile threshold is multiplied by 3 before being used for detection. Increasing this value makes ON detection more stringent; decreasing it increases sensitivity.
 - a. For now i have tested this parameter to be optimized for best detection (cf "Reglage du modèle")
11. The parameter `multiplier_tresh_accel` is a multiplicative factor applied to the acceleration threshold. For example, a value of 1 means the percentile is used directly. Values above 1 make Wake/REM detection more restrictive, while values below 1 increase sensitivity to movement.
 - a. For now i have tested this parameter to be optimized for best detection (cf "Reglage du modèle")
12. The parameter `saving` determines whether IS and Wake/REM timings are written to disk. When set to True, detected state transitions are stored at the end of the acquisition. When set to False, no file output is produced. This parameter has no influence on detection logic, only on data persistence.
 - a. DEPENDS ON YOU, if you want to save the results of the detection.
 - b. If True, be sure to indicate the path where it could be saved (cf 13.)
13. The parameter `path` defines the directory where output files are saved. For example, `"/media/data-103/Nathan/Avalanche_project/aval_data/"` specifies the storage location. The directory must exist and be writable; otherwise, file saving will fail.
 - a. DEPENDS ON YOU, if you want to save the results (the timings of IS) of the detection.
 - b. SOMETHING GOOD TO KNOW : i have called the file created : `"IS_wake_timings.txt"` (ligne 413 of the code), feel free to rename it on the code (ti change it : cf 14.
 - c. To be sure it will work well, create the file in advance to avoid a bug (happens at the very beginning i don't know why.
14. The parameter `file_name` is the name you choose to the file where the data of timings of the IS and wake/REM will be registered

15. The parameter `mean_spike_rate_channel` defines which channel index is used to read the mean spike rate signal. For example, channel 32 indicates that spike rate values are extracted from that channel in the Open-Ephys signal chain. An incorrect channel selection will result in invalid detection behavior.
 - a. DEPENDS ON YOU, really important to know where you have put the channel (Output of mean spike rate plugin)
16. The parameters `accel_channel_1`, `accel_channel_2`, and `accel_channel_3` define the channel indices corresponding to the three axes of the accelerometer. For example, 33, 34, and 35 indicate the channels used to compute total acceleration magnitude. If these channels are incorrectly assigned, Wake/REM detection will not reflect actual movement.
 - a. DEPENDS ON YOU, really important to know where you have put the channels (position numero in the list of the channel in the channel map, you can change their place in the channel map but you have to change it the numero of their place in this code too as well !)

SAVINGS :

At the end of the acquisition, the timings of IS and Wake/REM periods are saved in a text file named `IS_wake_timings.txt` in the directory defined by the parameter `path`. The file contains the acquisition number, the list of start markers for IS periods, and the list of start markers for Wake/REM periods. (cf PARAMETERS Section).

LAST ADVICES :

- Before starting any acquisition, the selected channels and thresholds should be adjusted if the recording setup differs from the original configuration.
- The TTL panel should be tested to ensure that `ligne_threshold_integration` and `ligne_stop_IS_now` trigger the expected behaviors.
- The saving directory should also be verified to confirm that the program has permission to write the output text file.
- if you do a change in the signal chain, i have always quit (on hubel with ctrl + C) and reopen open-ephys to refresh it (i know that is a shit functionality but that is what has worked until there lol). Because you want to be sure that your changes are taken into account

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2. **EXT**  **DAC**: External input mode is engaged when the **INPUT SELECT** switch is moved to the **EXT** position and user supplied voltage waveforms are present at the **EXT BNC** input. If the user attempts to supply more than 5V to the **EXT** input, the circuit will clamp the input signal to 5V. Using **EXT** mode, the behavior of the circuit is:

- **AUX**  **CURR**: Source the current specified by:

$$\text{Load Current} = \text{GAIN Position} * \frac{\text{SIGNAL IN Voltage}}{5V}$$

- **AUX**  **CURR**: Generate the optical power specified by:

$$\text{Power} = \text{GAIN Position} * h * \frac{\text{SIGNAL IN Voltage}}{5V}$$

The intensity of the LED will be dependent on the auxiliary feedback gain which defines the transmittance

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```

    def __init__(self, processor, num_channels, sample_rate):
        """

```

A new processor is initialized whenever the plugin settings are updated

Parameters:

processor (object): Python Processor class object used for adding events from python.

num_channels (int): number of input channels in the selected stream.

sample_rate (float): sample rate of the selected stream

```
        """
```

```
        print("Num Channels: ", num_channels, " | Sample Rate: ", sample_rate)
```

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```
def __init__(self, processor, num_channels, sample_rate):  
    """
```

A new processor is initialized whenever the plugin settings are updated

Parameters:

processor (object): Python Processor class object used for adding events from python.

num_channels (int): number of input channels in the selected stream.

sample_rate (float): sample rate of the selected stream

```
    """
```

```
print("Num Channels: ", num_channels, " | Sample Rate: ", sample_rate)
```

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```
self.start_time_cut = -1
```

```
self.processor.add_python_event(self.ttl_IS, False)
```

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```
avalanche_now = (current_mean-self.threshold_value) > 0
```

transition OFF -> ON (presence avalanche)

```
if avalanche_now and not self.avalanche_state:
```

```
    self.avalanche_state = True
```

```
    print("AVALANCHE ON | mean =", current_mean)
```

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```
if self.start_time_cut >= 0:
```

```
    self.end_time_cut = self.time_counter
```

```
    time_cut = self.end_time_cut - self.start_time_cut
```

```
    self.L_time_cuts.append(time_cut)
```

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```
    self.start_time_cut = -1
```

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```
self.processor.add_python_event(self.ttl_aval, True)
```

transition ON -> OFF (fin avalanche)

```
elif not avalanche_now and self.avalanche_state:
```

```
    self.avalanche_state = False
```

```
    print("AVALANCHE OFF | mean =", current_mean)
```

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```
if self.start_time_aval >= 0:
```

```
    self.end_time_aval = self.time_counter
```

```
    time_aval = self.end_time_aval - self.start_time_aval
```

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```

    def __init__(self, processor, num_channels, sample_rate):

```

```

        """

```

A new processor is initialized whenever the plugin settings are updated

Parameters:

processor (object): Python Processor class object used for adding events from python.

num_channels (int): number of input channels in the selected stream.

sample_rate (float): sample rate of the selected stream

```

        """

```

```

        print("Num Channels: ", num_channels, " | Sample Rate: ", sample_rate)

```

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```

        self.threshold_acceleration = None

```

```

        self.mean_acceleration_decay = None

```

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```

        self.mean_buffers = deque() # stocke (time, value)

```

```

        self.mean_buffer_time = 0

```

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```

        self.acceleration_buffers = deque() #stocke (time,value)

```

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pass

```

def handle_ttl_event(self, source_node, channel, sample_number, line, state):
    """
    Handle each incoming ttl event.

    Parameters:
    source_node (int): id of the processor this event was generated from
    channel (str): name of the event channel
    sample_number (int): sample number of the event
    line (int): the line on which event was generated (0-255)
    state (bool): event state True (ON) or False (OFF)
    """

    """# ici de 0 à x alors que sur open-ephys : de 1 à x+1
    if line == 24 and state == True: #marche pour 3,4,5,6,7,8 ... un peu près tous sauf 31
par (donc bien tester)
        self.processor.add_python_event(7, True) # marche pas pour 0,1,2 | marche pour
3,4,5,6,7
    if line == 24 and state == False:
        self.processor.add_python_event(7, False)"""

    ## CREATION D'UN TABLEAU DE CONTRÔLE avec TTL TOGGLE PANEL :
    # possibilite de lancer quand on le souhaite l'integration des threshold :
    if (line == self.ligne_threshold_integration and state == True) and not
self.user_redo_tresh:
        self.user_redo_tresh = True    # calcul du threshold demande
        self.user_counter_decision += 1 # +1 decision prise par le user sur le treshhold calcul
        self.processor.add_python_event(3, True)

    if (line == self.ligne_threshold_integration and state == False):
        self.processor.add_python_event(3, False)

    if line == self.ligne_stop_IS_now and state == True:
        self.user_stop_IS_now = True
        #self.processor.add_python_event(3, True)

```

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```
## transition ON -> OFF (fin avalanche)
elif not avalanche_now and self.avalanche_state:
    self.avalanche_state = False
    print("AVALANCHE OFF | mean =", current_mean)
```

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